

2. Modules and subworkflows

Nextflow from zero to hero

Modularity

and why it is your friend

- **Problem:** monolithic pipelines are hard to scale and maintain
- **Solution:** modularity
 - Maintainability
 - Code reuse across pipelines
 - Independent testing of components
 - Clear code organization
 - Community collaboration

Variant calling pipeline

FASTP

FASTQC

BWA_MEM

SAMTOOLS_SORT

SAMTOOLS_INDEX

SAMTOOLS_STATS

MOSDEPTH

DEEPVARIANT

MULTIQC

Scope of the “box”

Modules and Subworkflows

– Module

Process that is stored in a **separate file** and can be imported into different workflows.

– Subworkflow

Higher-level building block that groups **multiple modules** into a **logical unit of execution**, also in a separate file.

Variant calling pipeline

FASTP	]	READS_QC
FASTQC		
BWA_MEM	]	MAPPING
SAMTOOLS_SORT		
SAMTOOLS_INDEX		
SAMTOOLS_STATS	]	MAPPING_STATS
MOSDEPTH		
DEEPPVARIANT	]	VARIANT_CALL
MULTIQC	]	REPORTING

Creating a module

- Put the process in a separate, independent file with extension *.nf* (ie. *bwa.nf*)
- Input, logic and output of the module are controlled with `input`, `script` and `output` in the process definition.

```
process BWA_MEM {  
  
  input:  
  tuple val(sample_id), val(fastq_set_id), path(fastq_R1), path(fastq_R2), path(genome_fasta)  
  
  output:  
  tuple val(sample_id), path("${sample_id}-${fastq_set_id}.bwa.bam"), emit: bam_file  
  tuple val(sample_id), path("${sample_id}-${fastq_set_id}.bwa.log"), emit: bwa_log  
  
  script:  
  """"  
  bwa mem -t 4 \  
    -R "@RG\tID:${sample_id}\tSM:${sample_id}\tPL:Illumina" \  
    ${genome_fasta} \  
    ${fastq_R1} ${fastq_R2} \  
  2> ${sample_id}-${fastq_set_id}.bwa.log \  
    | samtools view --threads 4 -Sb - > ${sample_id}-${fastq_set_id}.bwa.bam  
  """"  
}
```

Creating a subworkflow

- **Import** the modules using `include`.
- Create a **named workflow**.
- Input, logic and output of the subworkflow are controlled with `take`, `main` and `emit` in the workflow definition.

```
include { BWA_MEM } from './bwa'  
include { SAMTOOLS_SORT } from './samtools-sort'  
include { MOSDEPTH } from './mosdepth'
```

Preview DAG

```
workflow MAPPING {  
  
    take:  
        qc_reads  
        genome_indexes  
  
    main:  
        bwa_input_ch = qc_reads.combine(genome_indexes)  
        BWA_MEM(bwa_input_ch)  
        bam_files_ch = BWA_MEM.out.bam_file  
  
        SAMTOOLS_SORT(bam_files_ch)  
  
    emit:  
        sorted_bam = SAMTOOLS_SORT.out.sorted_bam  
  
}
```